

# CleanTrust®

Fast. Accurate. Assured.



## What is ATP Testing?

ATP testing provides rapid, objective verification of hygiene after cleaning. By measuring residual organic material on surfaces or in water, ATP gives QA teams immediate feedback before production restarts. Results can be trended to show continuous improvement, support a proactive hygiene culture, and provide clear evidence during audits.

## ATP Surface Swabs

For verifying hygiene on food-contact and non-food-contact surfaces

### ATP Surface Swab

#### Key Features

- High sensitivity – detects ATP down to  $1 \times 10^{-15}$  mol (AOAC validated)
- High specificity
- Robust chemistry – minimal sanitiser interference
- Consistent performance – low lot-to-lot variation
- Operator-friendly
- Consistent RLU performance across multiple luminometers
- Compatible with CleanTrust Touch and Hygiena luminometers

#### Typical Applications

- Any food-contact surface
- Production equipment
- Conveyor belts
- Prep benches
- High-touch surfaces

#### Storage and Validity

- 12 months – store at 2–8°C
- 4 weeks – store at 20–25°C



Sanitiser-Resistant Performance

Compatible with CleanTrust Touch™ and Hygiena Luminometers

High-sensitivity ATP detection for critical hygiene monitoring



## ATP Surface MagniSwab

#### Key Features

- High sensitivity  $1 \times 10^{-16}$  mol
- Optimised swab design ensures maximum sample collection and recovery
- High-sensitivity, liquid-stable reagent enables detection of ultra-low ATP levels
- Excellent reproducibility
- Resistant to common sanitisers and disinfectants
- Fast and reliable results for real-time verification



# Water ATP Swabs

Total and Free ATP testing for CIP and water system hygiene.

## Key Features

- Sensitivity:  $1 \times 10^{-16}$  mol
- Total and Free ATP formats included
- Ideal for CIP rinse water and hard to reach surfaces
- Supports biocide checks, biofilm troubleshooting and water hygiene trending

## Storage and Validity

- 12 months – store at 2–8°C
- 4 weeks – store at 20–25°C



# CleanTrust Touch™ Luminometer

Backing Your Reputation with Irrefutable Cleaning Validation



## Key Features

- Touchscreen interface
- Stable optics for consistent RLU readings
- Calibration using purchased calibration rods
- Cloud-enabled data storage and reporting
- Results in seconds
- Rugged design for plant-floor use



## Workflow

Swab

Activate

Insert

Read

Record

## Ordering Information:

PRODUCT CODE	QTY	DESCRIPTION
AE0002	100	CleanTrust® ATP Surface Swab - Box
AE0003	100	CleanTrust® TOTAL ATP Water Swab - Box
AE0004	100	CleanTrust® FREE ATP Water Swab - Box
AE0012	100	CleanTrust® ATP MagniSwab - Box
AE1000	1	CleanTrust Touch™ Luminometer

1300 936 044

info@anaeron.com.au

20/8 Tilley Lane Frenchs Forest NSW 2086

www.anaeron.com.au

Fast. Accurate. Assured.

**ANAERON**  
Pty Ltd