




Simple Color Test

(No Luminometer Needed)



Lateral Flow Tests

(No Luminometer Needed)

Protein Residue Test

SpotCheck Plus

Glucose & Lactose

🕒 1 min

AllerSnap™

High-Sensitivity Broad Spectrum Protein

🕒 5-15 min

PRO-Clean

Protein

🕒 1-10 min

AlerTox® Sticks

Trace allergen contamination in food, liquid, and surfaces

Allergens that can be detected with AlerTox Sticks:



Crustacean



Egg



Hazelnut



Fish



β-Lactoglobulin



Almond



Mustard Seed



Casein



Soy



Total Milk



Walnut



Peanut

🕒 10 min

GlutenTox® Pro

AOAC Certified Gluten Detection harnessing the power of the G12 Antibody



🕒 10 min

GlutenTox® Sticks Plus

Detection limit 3 ppm (<20 ppm is considered "gluten-free")

🕒 10 min

Allerflow Gluten

High-Sensitivity
Gluten Specific

🕒 10 min

GlutenTox® Sticks Plus for Reader

Readable and trackable Quantitative or Pass/Fail results on the Hygiena Cube

🕒 15 min